User Manual

Wireless Arm Blood Pressure Monitor

Model: TMB-1585-BT

Distributed by

1byone Products Inc.

1230 E Belmont Street, Ontario, CA 91761

Customer Service: +1 909-391-3888

www.1byone.com

CE 0123

Guangdong Transtek Medical Electronics Co., Ltd.

Zone A, No.105, Dongli Road, Torch Development District,
Zhongshan, 528437, Guangdong, China

EC REP

MDSS - Medical Device Safety Service
GmbH Schiffgraben 41, 30175 Hannover, Germany

Thank you for selecting 1byone Wireless Arm Blood Pressure Monitor. Please read the instruction manual carefully and thoroughly so as to ensure the safe usage of this product. Keep this manual for further reference in case any issues arise.

Table of Contents

INTRODUCTION.....	2	DE
• General Description		
• Indication for Use		
• Safety Information		
• LCD Display		
• Product Description		
• Package List		
BEFORE YOU START.....	10	ES
• Choosing Power Supplies		
• Installing and Replacing the Batteries		
• Measurement Principles		
• Setting Date, Time and Measurement Unit		
• Selecting Users		
• Pairing the Blood Pressure Monitor with a Mobile Device		
MEASUREMENT.....	17	
• Fitting the Cuff		
• Starting the Measurement		
DATA MEASUREMENT.....	20	
• Viewing User Records		
• Deleting User Records		
GENERAL USER INFORMATION.....	23	
• Tips for Measurement		
• Maintenance		
ABOUT BLOOD PRESSURE.....	25	
• What is systolic pressure and diastolic pressure?		
• What is the standard blood pressure classification?		
• What is the irregular heartbeat detector?		
• Why does my blood pressure fluctuate throughout the day?		
• Why is my blood pressure different at home compared to the hospital?		
• Is the result the same if measuring on the right arm?		
TROUBLESHOOTING.....	28	
SPECIFICATIONS.....	29	
AUTHORISED COMPONENTS.....	30	
CONTACT INFORMATION.....	30	
COMPLIED STANDARDS LIST.....	31	
EMC GUIDANCE.....	32	
DISPOSAL.....	36	

EN General Description

This manual contains important safety and care information, and provides step by step instructions for using this blood pressure monitor. Read the manual thoroughly before using this product.

FR Blood Pressure Monitor Features

- 100 x 68mm (3.9" x 2.7") Digital LCD display with white backlight
- Maximum storage of 60 records per user
- Up-to-date measuring-during-inflation technology

IT Indication for Use

The 1byone Wireless Arm Blood Pressure Monitor is a digital monitor intended for use in measuring blood pressure and heartbeat rate with arm circumstances ranging from about 22 cm to 32 cm (about 8½"-12½"). It is intended for indoor, adult use only.

EN Safety Information

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "TYPE BF APPLIED PARTS"
	Symbol for "COMPLIES WITH MDD 93/42/EEC REQUIREMENTS"		
	Symbol for "MANUFACTURER"		Symbol for "ENVIRONMENT PROTECTION - Electrical waste products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "SERIAL NUMBER"		Symbol for "AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY"
	Symbol for "DIRECT CURRENT"		Caution: These notes must be observed to prevent any damage to the device.
	Symbol for "MANUFACTURE DATE"		Symbol for "RECYCLE"
	The Green Dot is the license symbol of a European network of industry-funded systems for recycling the packaging materials of consumer goods.		Symbol for "INCLUDING RF TRANSMITTER"

EN

CAUTION

- This device is intended for adult use in homes only.
- This device is not suitable for use on neonatal patients, pregnant women, patients with implanted, electronic devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral, arterial disease and patients undergoing intravascular therapy or arteriovenous shunt or people who received a mastectomy. Please consult your doctor prior to using this device if you suffer from illnesses.
- This device is not suitable for measuring the blood pressure of children. Ask your doctor before using it on older children.
- This device is not intended for patient transport outside a healthcare facility.
- This device is not intended for public use.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- Do not confuse self-monitoring with self-diagnosis. This device allows you to monitor your blood pressure. Do not begin or end medical treatment without asking a physician for treatment advice.
- If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your physician.
- Do not take any therapeutic measures on the basis of a self measurement. Never alter the dose of a medicine prescribed by a doctor. Consult your physician if you have any question about your blood pressure.
- When this device is used to measure patients who have common arrhythmia such as atrial or ventricular premature beats or atrial fibrillation, the best result may occur with deviation. Please consult your physician about the result.
- Don't kink or completely bend the connection tube during use, otherwise the cuff pressure may continuously increase which can prevent blood flow and result in harmful injury to the patient.
- When using this device, please pay attention to the following situations which may interrupt blood flow and influence blood circulation of the patient, thus causing harmful injury to the patient: kinking/bending the connection tube too frequently over consecutive measurements, applying and pressurizing the cuff on any arm undergoing intravascular access, therapy, or having an arteriovenous shunt present, or inflating the cuff on the same body side as a mastectomy.
- Warning: Do not apply the cuff over a wound; otherwise this may cause further injury.
- Do not inflate the cuff on the same limb which other monitoring ME equipment is applied around simultaneously, as this can cause temporary loss of function of the simultaneously-used monitoring ME equipment.

EN

CAUTION

- On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure (cuff pressure >300mmHg or constant pressure >15mmHg for more than 3 minutes) applied to the arm may lead to an ecchymosis.
- Please check that the operation of this device does not result in prolonged impairment of patient blood circulation.
- When measuring, please avoid compression or restriction of the connection tubing.
- This device cannot be used with HF surgical equipment at the same time.
- The accompanying documentation discloses the SPHYGMOMANOMETER was clinically investigated according to the requirements of ISO 81060-2:2013.
- To verify the calibration of the AUTOMATED SPHYGMOMANOMETER, please contact the manufacturer.
- This device is contraindicated for any female who may be suspected of, or is pregnant. Besides providing inaccurate readings, the effects of this device on the fetus are unknown.
- Taking measurements too frequently and consecutively with little time in between may cause disturbances in blood circulation and cause injuries.
- This device is not suitable for continuous monitoring during medical emergencies or operations. If used continuously, a patient's arm and fingers will become anaesthetic, swollen and even purple due to a lack of blood.
- When not in use, store this device and its adapter in a dry room and protect them against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on top of this device or its storage case.
- This device may only be used for purposes described in this manual. The manufacturer cannot be held liable for damage caused by incorrect application.
- This device is comprised of sensitive components and must be treated with caution. Observe the storage and operating conditions described in this manual.
- This device is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture of air with oxygen or nitrous oxide.
- Warning: No servicing/maintenance while the ME equipment is in use.
- The patient is an intended operator.
- The patient can measure, transmit data and change batteries under normal circumstances and maintain this device and its accessories according to the user manual.
- To avoid measurement errors, please avoid strong electromagnetic fields or fast transient/burst electrical signals as these can cause interference.
- This blood pressure monitor, the adapter and the cuff are suitable for use within the patient environment. If the patient/user is allergic to polyester, nylon or plastic, please don't use this device.

EN

CAUTION

- During use, the patient will be in contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.
- Adapter is specified as a part of ME EQUIPMENT.
- If you experience discomfort during measurement, such as pain in the arm or other complaints, press the START/STOP button to immediately release the air from the cuff. Loosen the cuff and remove it from your arm.
- If the cuff pressure reaches 40kPa (300mmHg), this device will automatically deflate. Should the cuff not deflate when pressures reaches 40kPa (300mmHg), detach the cuff from the arm and press the START/STOP button to stop inflation.
- Before use, make sure this device functions safely and is in proper working condition. Check this device. Do not use this device if it is damaged in any way. The continuous use of a damaged device may cause injury, improper results, or serious danger.
- Do not wash the cuff in a washing machine or dishwasher.
- The service life of the cuff may vary by the frequency of washing, skin condition, and storage state. Typical service life is 10,000 uses.
- It is recommended that performance should be checked every 2 years and after maintenance and repair by retesting at least the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50mmHg and 200mmHg).
- Please dispose of ACCESSORIES, detachable parts, and the ME EQUIPMENT according to the local guidelines.
- Manufacturer will make available on request circuit diagrams, component lists, descriptions, calibration instructions,etc., to assist service personnel in parts repair.
- The plug/adapter plug pins insulate this device from the power supply. Do not position this device in a position where it is difficult to disconnect from the power outlet to safely terminate operation of ME equipment.
- The operator shall not touch batteries /adapter and the patient simultaneously.
- Cleaning: Dusty environments may affect the performance of this device. Please use a soft cloth to clean the whole device before and after use. Don't use any abrasive or volatile cleaners.
- This device doesn't need to be calibrated within two years for reliable service.
- If you have any problems with this device, such as setting up, maintaining or using, please contact 1byone Customer Service. Don't open or repair this device by yourself in the event of malfunctions. This device must only be serviced, repaired and opened by individuals at authorized sales/service centers.
- Please report to 1byone if any unexpected operations or events occur.
- Keep this device out of reach of infants, young children or pets to avoid inhalation or swallowing of small parts. It is dangerous and may even be fatal.

EN

CAUTION

- Be careful of strangulation due to the length of this devices cables and hoses.
- At least 30 minutes is required for ME equipment to warm from the minimum storage temperature and be ready for use. At least 30 minutes is required for ME equipment to cool from the maximum storage temperature and be ready for use.
- This should be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, and walkietalkies can affect this device and should be kept at least a distance "d" away from the equipment. The distance "d" is calculated by the MANUFACTURER from the 800 MHz to 2.5 GHz column of Table 6 of IEC 60601-1-2:2007, as appropriate.
- Please only use ACCESSORIES and detachable parts specified/authorized by MANUFACTURER. Otherwise, unauthorized parts may cause damage to this device or danger to the user/patients.
- There are no Luer lock connectors used in the construction of tubing. There is a possibility that they might be inadvertently connected to intravascular fluid systems, allowing air to be pumped into a blood vessel.
- Please use this device in the appropriate environment as described in this user manual. Otherwise, the performance and lifetime of this device will be impacted and reduced.

DE

FR

ES

IT

DE

FR

ES

IT

EN ♥ LCD Display

DE

FR


ES

IT

Symbol	Description	Explanation
SYS.	Systolic blood pressure	High blood pressure
DIA.	Diastolic blood pressure	Low blood pressure
PUL/min	Pulse display	Pulse in beats per minute
▽	Deflation symbol	The cuff is deflating
M	Memory	Indicates it is in the memory mode
kPa	kPa	Measurement Unit of the blood pressure (1kPa=7.5mmHg)
mmHg	mmHg	Measurement Unit of the blood pressure (1mmHg=0.133kPa)
Lo +	Low battery	Batteries are low and need to be replaced
Heartbeat icon	Irregular heartbeat	Blood pressure monitor is detecting an irregular heartbeat during measurement
Bar chart icon	Blood pressure level indicator	Indicates the blood pressure level
PM MD 88/88	Current time	Year/Month/Day, Hour : Minute
Heartbeat icon	Heartbeat	Blood pressure monitor is detecting a heartbeat during measurement.
User A icon	User A	Start measurement and save and transmit the measuring results for User A
User B icon	User B	Start measurement and save and transmit the measuring results for User B
*	Bluetooth icon	The Bluetooth icon blinks when Bluetooth is working

FR ♥ Product Description

ES ♥ Package List

- Wireless Arm Pressure Monitor
- Cuff (Type BF applied part)

22 cm to 32 cm (about 8 3/4"-12 1/2")
Please only use the 1byone authorized cuff. Cuff sizing can be found on the cuff's label.

- Storage Bag


- Instruction Manual

EN ♥ Choosing Power Supplies

- DE** 1. Battery power mode:
6V DC 4 x AAA batteries (not included)
2. AC adapter power mode:
6V ---1A

FR Please only use the recommended AC adapter model. (not included)

ES Please unplug the adapter to disconnect from the using utility power.

IT CAUTION

In order to get the best from and protect your Arm Blood Pressure Monitor, please use the correct batteries and power adapters which comply with CE safety standards.

♥ Installing and Replacing the Batteries

1. Open the battery cover.
2. Install the batteries by matching the correct polarity.
3. Reinsert the battery cover.

Replace the batteries when:

- The **Lo +** indicator is shown.
- The display dims
- The display does not light up

IT CAUTION

- Do not use new and used batteries together.
- Do not use different types of batteries together.
- Remove batteries if the device is not likely to be used for a long period of time.
- Old batteries are harmful to the environment. Do not dispose with other daily trash.
- Remove the old batteries from the device and follow your local recycling guidelines.
- Do not dispose of batteries in fire. Batteries may explode or leak.

EN ♥ Measurement Principles

This device uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the device establishes a "zero pressure" equivalent to the air pressure. Then it starts inflating the arm cuff. Meanwhile, the device detects pressure oscillations generated by beat-to-beat pulsatile, which are used to determine the systolic and diastolic pressure, and also pulse rate.

The device also compares the longest and the shortest time intervals of detected pulse waves to the mean time interval, then calculates the standard deviation. The device will display a warning signal with the reading to indicate the detection of an irregular heartbeat when the difference of the time intervals is over 25%.

♥ Setting Date, Time and Measurement Unit

It is important to set the clock before using your Arm Blood Pressure Monitor, so that a time stamp can be assigned to each record that is stored in the memory. Year range: 2014-2054. Time formats: 12H/24H.

1. When the Arm Blood Pressure Monitor is off, press and hold the "S" button for 3 seconds to enter the year setting mode.

2. Press the "M" button to change the [YEAR]. Each press will increase the number by one in a cycling manner. Press and hold the "M" button to increase the number faster.

BEFORE YOU START

EN


3. Press the "S" button to confirm the year and go to the next step. You can switch between the [DATE FORMAT] of D/M or M/D by pressing the M button. Press and hold the "M" button to increase the number faster.

DE


FR

ES


IT

4. Repeat steps 2 and 3 to set the [MONTH] and [DAY].


5. Repeat steps 2 and 3 to set the [TIME FORMAT] between 12H and 24H.


BEFORE YOU START

EN

6. Repeat steps 2 and 3 to set the [HOUR] and [MINUTE].

7. Repeat steps 2 and 3 to set the [UNIT].

8. After the device is set, the LCD will display "dONE" and all the settings will be displayed. Then the display will turn off.


EN

♥ Selecting Users

DE

- When the Arm Blood Pressure Monitor is off, press and hold the "M" button for 2 seconds to enter user setting mode. The user ID will blink.


FR

ES

IT

- Press the "M" button again to select the user ID between user A and user B.

- After selecting the suitable user ID, press the "S" button to confirm. The display will then turn off.


♥ Pairing the Blood Pressure Monitor with a Mobile Device

EN


DE

Before pairing, please download the MedM Health App from Apple App Store or Google Play Store.

- Turn on Bluetooth and the App. Make sure both are ON.
- When the monitor is OFF, press and hold the "START/STOP" button to start pairing. The Bluetooth symbol * will blink, indicating pairing is proceeding.

- Tap the icon on the MedM Health app and choose "My Devices" to add a new device.

EN 4. Select the Smart BP Monitor 1585B and then select the user ID you want to connect with. Select User profile 1 for User A and User profile 2 for User B.

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If successful, symbol will be shown on the LCD.

If pairing has failed, only Bluetooth symbol will be shown on the LCD.

5. The monitor will shut off after the pairing process is complete.

Note: After pairing is complete, you can start a measurement and data will synchronize with the App automatically in a few seconds.

Bluetooth Module No.: AW2540MV1
RF Frequency Range: 2402 MHz to 2480 MHz
Output Power Range: ≤-1 dBm
Supply Voltage: 2V-3.6 V
Transmitting Distance: 10 meters

List of compatible devices

- For iOS devices:
The operating system must be iOS 8 or higher, as available on iPhone 4S, iPhone 5/5C/5S, iPhone 6/6 Plus and so on.
- For Android devices:
The operating system must be 4.3 or higher.

Fit the Cuff

1. Remove all jewelry, such as watches and bracelets from your left arm.

Note: If your doctor has diagnosed you with poor circulation in your left arm, use your right arm.

2. Roll or push up your sleeve to expose the skin. Make sure your sleeve is not too tight.

3. Hold your arm with your palm facing up and tie the cuff on your upper arm, then position the tube off-center toward the inner side of arm in line with the little finger. Position the artery mark over the main artery (on the inside of your arm).

Note: Locate the main artery by pressing with 2 fingers approximately 0.79in (2cm) above the bend of your elbow on the inside of your left arm. Identify where the pulse can be felt the strongest. This is your main artery.

4. The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.

5. Sit comfortably with the monitored arm resting on a flat surface. Place your elbow on a table so that the cuff is at the same level as your heart. Turn your palm upwards. Sit upright in a chair, and take 5-6 deep breaths.

6. Helpful tips for Patients, especially for Patients with Hypertension:

- Rest for 5 minutes before first measuring.
- Wait at least 3 minutes between measuring. This allows your blood circulation to recover.
- Take the measurement in a silent room.
- The patient must relax as much as possible and should not move and talk during the measurement procedure.
- The cuff should maintain its position at the same level as the right atrium of the heart.
- Please sit comfortably. Do not cross your legs and keep your feet flat on the ground.
- Keep your back against the backrest of the chair.
- For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.


EN ♥ Starting the Measurement

DE 1. When the monitor is off, press the "START/STOP" button to turn on the monitor and complete a full measurement. (Below example is for User A)


FR LCD display


ES Inflating and measuring

IT Adjusting to zero

DE Final, saved measurement results

FR 2. This device will process and transmit after measurement. The Bluetooth symbol will blink on the LCD indicating data is transmitting.

ES If data transmission fails, the LCD will display **Err**.

IT If data transmission succeeds, the LCD will display **d0 RE**.

EN 3. Press the "START/STOP" button to power off.

Tips: Maximum 60 records each will be stored for User A and User B.

CAUTION

- Interference may occur in the vicinity of equipment marked with the following symbol: . TMB-1585-BT interference may also occur in the vicinity of electrical equipment.
- Sensitive people, including pregnant women pre-eclamptic and those with implanted medical electronic instruments, should avoid using the device whenever possible.
- Keep the monitor at least 20 centimeters away from the human body (especially the head) when data transmission is proceeding after measurement.
- To enable the data transmission function, this product should be paired to Bluetooth of at least 2.4 GHz.

How to mitigate possible interference?

1. The range between the monitor and mobile BT device should be reasonably close, from 1 meter to 10 meters. Please ensure no obstacles between the monitor and mobile BT device so as to obtain quality connection and to lower the RF output range.
2. Avoid interference, other electronic devices (particularly those with wireless transmission or a transmitter) should be kept at least 1 meter away from the monitor.

EN Viewing User Records

DE 1. When the Arm Blood Pressure Monitor is off, please press the "M" button to show the average value of the latest three records.

FR Note: If less than 3 records are stored, the LCD will display the most recent record instead. See User A on the right as an example.

ES

IT 2. Press the "M" button or the "S" button to view the record you want. The order, date and time of the record will be shown alternately.

3. If you want to check another user's records, press the "START/STOP" button to turn off the Arm Blood Pressure Monitor when it is in memory inquiry mode. Press and hold the "M" button for 2 seconds to enter into the selecting user ID mode, press the "M" button again to select the user ID between user A and user B, press the "S" button to confirm the user ID and then press the "M" button to check the selected user's measurement records.



CAUTION

The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (60) is dropped from the list.

EN 4. Users can also view records on the App in List, view by day, week or month.

August 17, 2017	
16:05	88 / 57 (53)
11:40	96 / 58 (60)
11:09	94 / 67 (61)
09:15	119 / 50 (96)

EN ♥ Deleting User Records

If there has been a measurement error or a correct measurement is not taken, delete the result by following the steps below. (Example shown below for User A)

- DE** 1. Press and hold the "M" button for 3 seconds when the Arm Blood Pressure Monitor is in the memory recall mode. The display will flash "User ID+del ALL."
- FR**
- ES**
- IT**


2. Press the "S" button to confirm deletion and the Arm Blood Pressure Monitor will display "dEL dOnE" and then turn off.

Note: To exit the delete mode without deleting any records, press the "START/STOP" button before pressing the "S" button.

3. If there is no record, the following display will be shown.

EN ♥ Tips for Measurement

Measurements may be inaccurate if taken in the following circumstances.

EN

♥ Maintenance

In order to get the best performance, please follow the instructions below.

DE

Put in a dry place and avoid sunshine.

Avoid intense shaking and collisions.

Using wet cloths to remove dirt.

FR

Avoid touching water. Clean any water with a dry cloth.

Avoid dusty and unstable temperature environments.


Do not attempt to clean the reusable cuff with water and never immerse the cuff in water.

ES

IT

♥ What is systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.

EN

DE

FR

ES


IT

♥ What is the standard blood pressure classification?

The blood pressure classification published by the World Health Organization (WHO) and International Society of hypertension (ISH) in 1999 is as follows:

⚠ CAUTION -

Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

Blood Pressure (mmHg)	Level	Optimal	Normal	High-normal	Mild	Moderate	Severe
SYS	<120	120-129	130-139	140-159	160-179	≥180	
DIA	<80	80-84	85-89	90-99	100-109	≥110	

EN ♥ What is the irregular heartbeat detector?

DE An irregular heartbeat is detected when a heartbeat rhythm varies while the device is measuring the systolic and diastolic blood pressure. During each measurement, this equipment records the heartbeat intervals and works out the standard deviation. If the calculated value is larger than or equal to 25%, the irregular heartbeat symbol appears on the display when the measurement results are displayed.

FR

ES**IT**

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heart-beat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

♥ Why does my blood pressure fluctuate throughout the day?

1. Individual blood pressure varies multiple times everyday. It is also affected by the way the cuff is fitted and your measurement position, so please take the measurement under the same conditions.
2. Taking medication may cause measurement results to vary more widely.
3. Wait at least 3 minutes between measurements.

♥ Why is my blood pressure different at home compared to the hospital?

EN Pay attention to the following when measuring blood pressure at home:


Blood pressure may be different even throughout the day due to weather, emotion, exercise etc. Also, there is the "white coat" effect, which means blood pressure usually increases in clinical settings.

- The cuff is properly fitted.
- The cuff isn't too tight or loose.
- The cuff is fitted to the upper arm.
- User is relaxed and does not feel anxiety.
- User relaxes for 4-5 minutes before measurement.
- User takes 2-3 breaths before measurement begins.

♥ Is the result the same if measuring on the right arm?

IT

It is OK for both arms, but there will be some different results for different people. We suggest you measure the same arm every time.

TROUBLESHOOTING

SPECIFICATIONS

EN This section includes a list of error messages and frequently asked questions for problems you may encounter with your Arm Blood Pressure Monitor. If the product is not operating as you think it should, check here before arranging for servicing.


	Problem	Symptom	Check This	Remedy
FR ES IT	No Power	Display will not light up.	Batteries are exhausted.	Replace with new batteries.
			Batteries are inserted incorrectly.	Insert the batteries correctly.
	Low Batteries	Display is dim or show  +  .	Batteries are low.	Replace with new batteries.
	Error Message	 +Err shows.	Data communication/transmission failed.	Check that App/Bluetooth is on. Try data transmission again.
		E1 shows.	The cuff is not secure.	Refasten the cuff and then measure again.
		E2 shows.	The cuff is very tight.	Refasten the cuff and then measure again.
		E3 shows.	The cuff has excess pressure.	Relax for a moment and then measure again.
		E10 or E11 shows	The monitor detects motion or talking, or the user's pulse is too weak.	Relax for a moment and then measure again.
		E20 shows	The monitor does not detect the pulse signal.	Loosen the clothing on the arm and then measure again.
		E21 shows	The measurement failed.	Relax for a moment and then measure again.
		EExx shows on the display.	A calibration error occurred.	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions.

Power Supply	Battery powered mode: 6V DC 4 x AAA batteries (not included) AC adapter powered mode: 6V  1A Please use the recommended AC adapter model (not included)
Display Mode	Digital LCD with white backlight V.A. 100 x 68mm (3.9 x 2.7")
Measurement Mode	Oscillographic testing mode
Measurement Range	Rated cuff pressure: 0mmHg-299mmHg(0kPa ~ 39.9kPa) Measurement pressure: SYS: 60mmHg-230mmHg (8.0kPa~30.7kPa) DIA: 40mmHg-130mmHg (5.3kPa~17.3kPa) Pulse value: (40-199)beat/minute
Accuracy	Pressure: 41°F-104°F (5°C-40°C) within±3mmHg (0.4kPa) Pulse value: ±5%
Normal working condition	Temperature: 41°F to 104°F (5°C to 40°C) Relative Humidity: ≤85%RH Atmospheric pressure: 86kPa to 106kPa
Storage & transportation condition	Temperature: -4°F to 140°F (-20°C to 60°C) Relative Humidity: 10%RH-93%RH Atmospheric Pressure: 50kPa-106 kPa
Measurement perimeter of the upper arm	About 22 cm to 32 cm (about 8½"-12½")
Net Weight	Approx.282g (Excluding the dry cells and cuff)
External dimensions	Approx.154 x 106 x 57.1mm (6.1 x 4.2 x 2.2")
Attachment	1 x storage bag, instruction manual
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Protection against ingress of water	IP21 means monitor is protected against solid foreign objects of 12.5mm and greater, and protect against vertically falling water drops.
Device Classification	Battery Powered Mode: Internally Powered ME Equipment AC Adapter Powered Mode: Class II ME Equipment
Software Version	V01

WARNING: No modification of this equipment is allowed.

AUTHORIZED COMPONENTS

COMPLIED STANDARDS LIST

2. Storage Bag

Risk management	EN ISO 14971:2012 / ISO 14971:2007 Medical devices Application of risk management to medical devices
Labeling	EN 980:2008 Symbols for use in the labelling of medical devices ISO 15223-1:2012 Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 1: General requirements
User manual	EN 1041:2008 Information supplied by the manufacturer of medical devices
General Requirements for Safety	EN 60601-1:2006 / IEC 60601-1:2005+A1:2012 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance EN 60601-1-11:2010/ IEC 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Electromagnetic compatibility	EN 60601-1-2:2007/ IEC 60601-1-2:2007 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
Performance requirements	EN ISO 81060-1:2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type EN 1060-3:1997+A2:2009 Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems IEC 80601-2-30:2013 Medical electrical equipment- Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers
Clinical investigation	EN 1060-4:2004 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers ISO 81060-2:2013 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated measurement type
Usability	EN 60601-1-6:2010/IEC 60601-1-6:2010+A1:2013 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability EN 62366:2008/ IEC 62366-1:2015 Medical devices - Application of usability engineering to medical devices
Software life-cycle processes	EN 62304:2006/AC: 2008 / IEC 62304:2006 Medical device software - Software life-cycle processes
Bio-compatibility	ISO 10993-1:2009 Biological evaluation of medical devices- Part 1: Evaluation and testing within a risk management process ISO 10993-5:2009 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity ISO 10993-10:2010 Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization

EN EMC Statement

- DE 1) This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this device can be affected by portable and mobile RF communications equipment.
- FR 2) Do not use a mobile phone or other devices that emit electromagnetic fields, near the device. This may result in incorrect operation of the device.
- ES 3) Caution: This device has been thoroughly tested and inspected to assure proper performance and operation!
- IT 4) Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Table 1 Guidance and MANUFACTURER's declaration – ELECTROMAGNETIC EMISSIONS- for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emissions		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2 Guidance and MANUFACTURER's declaration – electromagnetic IMMUNITY – for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

EN Table 4 Guidance and MANUFACTURER's declaration – electromagnetic IMMUNITY – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

DE **FR** **ES** **IT** **Guidance and manufacturer's declaration – electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}$ 80MHz to 800 MHz $d=2,3\sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

EN **DE** **FR** **ES** **IT** **Table 6 Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 25 GHz
$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

EN Disposal of the appliance

Under no circumstances should you dispose of the appliance in normal domestic waste. This product is subject to the provisions of European Directive 2012/19/EU.

FR Dispose of the appliance via an approved disposal company or your municipal waste facility. Please observe currently applicable regulations. Please contact your waste disposal center if you need any further information.

ES Disposal of batteries

Used batteries may not be disposed of in household waste.

All consumers are statutorily required to dispose of batteries at a collection point provided by their local municipality or retail store.

The purpose of this obligation is to ensure batteries are disposed of in a non-polluting manner. Only dispose of batteries when they are fully discharged.

The appliance's packaging is made from environmentally friendly material and can be disposed of at your local recycling plant.

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IT Kontakt Informationen

For more information about our products, please visit www.1byone.com. Users can also connect with customer service and receive customer downloads and support.

1byone Products Inc.
1230 E Belmont Street, Ontario, CA 91761
Customer Service:+1 909-391-3888
www.1byone.com

Manufactured by: Guangdong Transtek Medical Electronics Co., Ltd.
Company: Guangdong Transtek Medical Electronics Co., Ltd.
Address: Zone A, No.105, Dongli Road, Torch Development District, Zhongshan, 528437, Guangdong, China

Authorized European Representative:
Company: MDSS - Medical Device Safety Service
Address: GmbH Schiffgraben 41, 30175 Hannover, Germany

Risk management	EN ISO 14971:2012 / ISO 14971:2007 Medical devices Application of risk management to medical devices
Labeling	EN 980:2008 Symbols for use in the labelling of medical devices ISO 15223-1:2012 Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 1: General requirements
User manual	EN 1041:2008 Information supplied by the manufacturer of medical devices
General Requirements for Safety	EN 60601-1:2006 / IEC 60601-1:2005+A1:2012 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance EN 60601-1-11:2010 / IEC 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Electromagnetic compatibility	EN 60601-1-2:2007 / IEC 60601-1-2:2007 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
Performance requirements	EN ISO 81060-1:2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type EN 1060-3:1997+A2:2009 Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems IEC 80601-2-30:2013 Medical electrical equipment- Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers
Clinical investigation	EN 1060-4:2004 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers ISO 81060-2:2013 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated measurement type
Usability	EN 60601-1-6:2010/IEC 60601-1-6:2010+A1:2013 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability EN 62366:2008 / IEC 62366-1:2015 Medical devices - Application of usability engineering to medical devices
Software life-cycle processes	EN 62304:2006/AC: 2008 / IEC 62304:2006 Medical device software - Software life-cycle processes
Bio-compatibility	ISO 10993-1:2009 Biological evaluation of medical devices- Part 1: Evaluation and testing within a risk management process ISO 10993-5:2009 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity ISO 10993-10:2010 Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization

EN EMC Statement

- 1) This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this device can be affected by portable and mobile RF communications equipment.
- 2) Do not use a mobile phone or other devices that emit electromagnetic fields, near the device. This may result in incorrect operation of the device.
- 3) Caution: This device has been thoroughly tested and inspected to assure proper performance and operation!
- 4) Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

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DE Guidance and manufacturer's declaration – electromagnetic immunity

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FR	Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
ES	Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
IT	Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Recommended separation distance $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}$ 80MHz to 800 MHz $d=2,3\sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.b Interference may occur in the vicinity of equipment marked with the following symbol: 

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EN Table 6 Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

DE Recommended separation distances between portable and mobile RF communications equipment and the device.

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FR Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 25 GHz
ES	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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