# SMART SCALE

1byone Products Inc.  
1230 E Belmont Street, Ontario, CA 91761  
Customer Service: +1 909-391-3888  
[www.1byone.com](http://www.1byone.com)

Instruction Manual/Bedienungsanleitung/Manuel d'instructions  
Manual de instrucciones/Manuale di istruzioni

Thank you for purchasing the 1byone Smart Scale. Please read this instruction carefully before operating.

## Measure Range

Body weight, Body fat, Body water, Muscle mass, BMI(Body mass index), BMR(KCAL), Bone mass and Visceral fat.

## Working Theory

This device is based on BIA (Bioelectrical Impedance Analysis) technology. To measure body bioelectrical impedance, it makes the results much more accurate with new measuring technology. The reason is that human body fat is not conductive, but muscle and water are good conductors. The device measures human body bioelectrical impedance to calculate the composition of human body.

## Caution

### Safety Precautions

1. Not recommended for any person using a pacemaker or other internal device.
2. Not recommended for women who are pregnant.
3. Any measurement obtained using this device is for reference only and should not be considered as a medical opinion.
4. Consult your primary care doctor or physician before making changes to your diet, exercise plans or physical activities.
5. DO NOT stand on the edge of the scale or jump on it.
6. DO NOT overload the scale (Maximum 396lb/180kg).
7. Always place scale on a hard, dry and flat surface before measurement.
8. Make sure your feet are dry before stepping on the scale.
9. Connect the batteries in accordance with the correct polarities indicated.
10. Remove batteries if the scale is not to be used for a long time.

**Note**

- Do not bend knees when measuring.
- Take measurements at the same or similar time of the day.
- Do not take measurements immediately after intense exercise.
- Mobile phone and microwave frequencies may affect accuracy of results.
- The result is just a reference for children, aged people and professional athletes.
- Do not disassemble it by yourself.
- Do not wash the scale under the water.
- Do not wash it with hot water or oil.
- Wiping the machine by wet cloth with detergent, and clean it with dry cloth.
- Before using the scale, make sure you have programmed all your personal data.
- When you enter "Baby Mode", please weigh with shoes.
- Take measurements at the same or similar time of the day.
- You can change the unit by "Setting>My information>Unit".

**Preparation Before Use****Inserting Batteries**

To insert the batteries, open the battery compartment on the rear of the scale and remove the compartment cover. Insert 4 AAA batteries, paying attention to the correct polarity.

**Auto Power-off**

The scale will automatically power off:

- After data has been uploaded to the application.
- After 30 seconds of no operation.
- It is overloaded.
- Due to faulty operations.

**Switch unit**

1. Three units: kg/lb/st:lb for option.
2. The default unit on scale is decided by the UNIT Selection in APP, please change default unit in APP for daily use.  
Please note that if user select "st:lb" in APP, the scale will show "lb" only.

**Operation Instruction****What devices are compatible?****iOS devices:**

Phone 4s, iPad3, iPad Mini or higher.

**Android devices:**

Android 4.3 or higher with Bluetooth 4.0 or higher version, memory capacity no less than 1GB.

**How to install the 1byone Wellness 2.0 APP**

Download App "1byone Wellness 2.0" from App Store or Google Play.

1. Make sure your phone or tablet is connected to the Internet.
2. Search and download "1byone Wellness 2.0" in App Store or Google Play. Or simply scan the QR code below.

iOS  
1byone Wellness 2.0

Android  
1byone Wellness 2.0

1byone Wellness 2.0  
APP icon

3. The APP icon will appear at the desktop of your phone or tablet after the installation is completed.
4. Tap "1byone Wellness 2.0" icon to launch the App.

## How to use the 1byone Wellness APP?

1. Turn on the Bluetooth function of your device.
2. Tap "1byone Wellness 2.0" to open the APP on your device.
3. When using the APP for the first time, you need to follow the steps below to enter your Profile information.

A. This step only applies to iPhone that support iOS 8 or higher. If this doesn't apply to you, please go to step B. When you open 1byone Wellness 2.0 App for the first time, the system will automatically turn to health data access page. You can select all health data types which can be synchronized to Health, and tap "Done" on top right corner afterwards. Currently the Health App only accepts 5 parameters: Body Fat Percentage, BMI, Height, Lean Body Mass and Weight.

B. As you are prompted, accurately fill in personal information. Tap on the picture to change user's image. Input your information (user name, gender, preferred weight unit, birthday, height, fitness level), check and confirm, then save to continue.

1

2

C. You will see the following indication, keep standing on the metal sheets of the scale with bare feet to make a connection between your device and scale via Bluetooth.

3

4

D. After the scale display blinks "---" three times, it will present the confirmed weight measurement. That means the body composition analysis are measured, click "Yes" to receive them to your app. You can check the standard of your each body composition by tapping on the data.

5

6

E. Click into the history, you can see your history data in graph or list. You can share your body data via Twitter/Email/SMS/Facebook.

7

Date	Weight	Stat
19-11-14, 09:00:00	10.10kg	10.1
19-11-14, 09:00:00	10.05kg	10.1
19-11-14, 09:00:00	10.05kg	10.1
19-11-14, 09:00:00	10.05kg	10.1

8

9

## Baby Mode

1. To create an account for you baby, click the Baby mode and input the baby's information and click "Save".

10

Name:   
 Gender:  Male  Female  
 Length:  cm  
 Birthday:   
 Weight:  kg  
 Name:  kg

11

06

2. Click "Weighing" and stand on the scale with your shoes when the scale showing "Adult weighing first", leave the scale and stand on it again with your baby when the indication turns to "Weighing with baby".

12

13

3. After weighing click "Yes" to save the data. You can also see the history data in graph and list by pressing history.

14

Date	Weight	Stat
19-11-14, 09:00:00	10.10kg	10.1
19-11-14, 09:00:00	10.05kg	10.1
19-11-14, 09:00:00	10.05kg	10.1
19-11-14, 09:00:00	10.05kg	10.1

15

07

## Navigating the Main Menu

16

1. Add and delete user: tap the user image icon  to enter the "User Select" page.
  - Tap "+" to add new users, each user needs to re-input user information.
  - Tap "Edit", then tap "-" to remove users.
2. History: measured data
3. Settings: Reset user information, save data in TXT files, information about body health reference tables and view the application's description.
4. Baby mode


## Compatible with Apple "Health" APP

There is a health APP called Health which is compatible to 1byone Wellness 2.0. You can update body data from 1byone Wellness 2.0 to Health. You can check historical data even if you uninstall 1byone Wellness 2.0 accidentally.

Notes: For iOS8 or higher.


## How to Use Health App?

1. Find Health App on your iPhone.
2. When using 1byone Wellness for the first time, you will be asked for permission to update date from 1byone Wellness 2.0 to Health.

17

3. Please note that date is only updated in Health when the measurements are successfully taken in 1byone Wellness 2.0. Unsuccessful measurements will not be updated to Health.
4. Open Health App, choose the "Body Measurements", you will see the latest data updated from 1byone Wellness 2.0.

18

## Frequently Asked Questions & Answers

- Why doesn't the scale work? Why doesn't the display show data? Why does the data on the screen disappear in a flash?
  - Please check whether the insulation piece is taken out, batteries are installed or batteries are in low and need replacing.
- Display indicates "Err" when measuring body fat.
  - Your feet or scale are too wet.
  - Overload the maximum capacity (396lb/180kg)
  - Body parameter settings should be with Height(3'3"-7'3"/100-220cm), Aged(10-99 years old).
- The scale can not connect with App.
  - The app is required to be opened when weighing.
  - Make sure to turn on the Bluetooth function on your device.
  - Please re-turn on the Bluetooth and restart 1byone Wellness 2.0 on your device.
  - Check the phone model, the phone must be higher than iPhone 4S, iPad 3, or iPad Mini for iOS. For Android models, phones must be Android 4.3 or higher.
- No body fat measured when weighing.
  - Body fat scale requires weighing with dry and bare feet.
  - Ensure the Bluetooth function is open.
  - Please input your personal information with Height(3'3"-7'3"/100-220cm), Aged(10-99 years old).

## Maintenance and Cleaning

- Do not keep the scale loaded for a long time. This may cause damage to the scale and it will no longer be accurate.
- Do not put anything on the scale when not in use. Otherwise, it will consume power and reduce the lifetime of the scale and batteries.

- Keep the product clean and dry; Keep away from water; Avoid damaging the inside electric circuit when weighing.
- Do not drop the product. Avoid risking damaging the scale or breaking the glass.
- The scale is in lower power when "Lo" is displayed. Please change the batteries immediately. Permanent damage will occur if the scale is kept in a low power mode for a long time.
- This product is a body analyzer scale for body composition analyzing only. Please do not use for commercial purposes.
- Clean with a damp cloth and wipe down using a dry cloth.
- DO NOT use detergents/cleaning agents to clean the scale. Doing so may permanently damage the plating on the electrodes used to measure body fat.

## Accuracy Range

50KG:±300g	100KG:±400g	150KG:±500g
110lbs+/-0.66lbs	220lbs+/-0.88lbs	330lbs+/-1.1lbs

## Indication Symbols

Lo	Low battery indicator
Err	Overload indicator

## Specifications

- 27x80mm Highlight LED Display
- Unit: kg/lb/st:lb
- 10 user groups
- Results are displayed graphically
- Three fitness levels: Ordinary, Amateur and Professional
- Capacity:180kg/396lb
- Minimum Mass measurement graduation: 0.1kg
- Muscle Mass measurement graduation: 0.1kg
- Weight measurement graduation: 0.1kg/0.2lb

- EN 10. Body Fat measurement graduation: 0.1%
- DE 11. Body Water measurement graduation:0.1%
- FR 12. Bone graduation: 0.1kg
- ES 13. Age range: 0~99 years old
- IT 14. Height range: 100~220cm 3ft:3in~7ft:3in
- 15. Working temperature:10~40°C
- 16. Working humidity:20%~90%
- 17. Storage temperature: -10~50°C
- 18. Power: DC 6V (4 x AAA Batteries)

### DISPOSAL

#### Disposal of the appliance

Under no circumstances should you dispose of the appliance in normal domestic waste. This product is subject to the provisions of European Directive 2012/19/EU.

Dispose of the appliance via an approved disposal company or your municipal waste facility. Please observe currently applicable regulations. Please contact your waste disposal center if you need any further information.

#### Disposal of batteries

Used batteries may not be disposed of in household waste.

All consumers are statutorily required to dispose of batteries at a collection point provided by their local municipality or retail store. The purpose of this obligation is to ensure batteries are disposed of in a non-polluting manner. Only dispose of batteries when they are fully discharged.

The appliance's packaging is made from environmentally friendly material and can be disposed of at your local recycling plant.

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